

- REF** 80538 INNO-LIA HCV Score
- REF** 80540 INNO-LIA HIV I/II Score
- REF** 80541 INNO-LIA HTLV I/II Score
- REF** 80542 INNO-LIA Syphilis Score

**Manufactured by:**

Fujirebio Europe N.V.
Technologiepark 6
9052 Gent
Belgium
☎ +32 9 329 13 29
BTW BE 0427.550.660
RPR Gent

Summarized Test Procedure in addition to the instructions for use of the INNO-LIA Score products

Note changes highlighted

Distributed by:

Fujirebio Europe N.V.
☎ +32 9 329 13 29
customer.support@fujirebio.com

Fujirebio Italia S.r.l.
☎ +39 06 965 28 700
italy@fujirebio.com

Fujirebio Germany GmbH
☎ +49 511 857 3931
germany@fujirebio.com

Fujirebio Iberia S.L.
☎ +34 93 270 53 00
spain@fujirebio.com

Fujirebio France SARL
☎ +33 1 69 07 48 34
france@fujirebio.com



eifu.fujirebio.com/FRE/all/?keycode=FRI00560


☎ +800 135 79 135 / +31 20 794 7071

8:00 – 17:00 GMT+1

M	T	W	T	F	S	S
☒	☒	☒	☒	☒	☐	☐

Summarized Test Procedure

ROBOBLOT




This summarized test procedure is a quick guide to perform an INNO-LIA Score test on the *RoboBlot*.

Before using the instrument for the first time, read the detailed “RoboBlot User Manual” which is provided with every instrument.


I. Setup of Instrument

1. Switch the instrument ON

- Use the ON/OFF switch on the right side of the instrument to switch the instrument ON
- Wait until the initialization is done and the instrument displays:



- Press Start





- The RoboBlot is in standby mode and ready to start an INNO-LIA test

II. Reagent Preparation and Start Run



1. Select a test procedure

- Bring the reagents to room temperature (see instructions for use of the assay)
- When starting a test manually on the RoboBlot, press the “Start Assay” key
- Type the name of your test in the “Worklist Name” field
- Select the assay protocol you wish to run



- Select the Starting well on the tray from where the assay will begin

- Select the Number of Samples and indicate if you want to run the Positive and/or Negative control





REMARK:


- The numbering on the tray starts at 0, but the starting well on the screen is the absolute position in the tray.
e.g.: Starting well 5 = number 4 on the tray

• When starting a Run that was prepared in the LiRAS software, select “Assays Worklists” – “Pre-Generated Worklists” and indicate your pre-generated worklist. Follow the instructions on screen as described below to start the run

- Configure the first tip picking position and fill the rack as shown. Use only dedicated tips (see references)




- The volumes required for this run are calculated by the instrument. Fill the dedicated bottles with minimum the amount of reagents indicated in the software



- Insufficient reagent will pause the instrument during the run

2. Prepare the reagents

- Prepare the Drying solution, a 20% (V/V) EtOH solution
- Fill the original, dedicated bottles with the correct reagents and place them on the corresponding spaces in the reagent holders, minimum volumes are indicated on screen
- Add the Positive and Negative control bottles on the dedicated place in the reagent holder



For INNO-LIA Score assays:

INNO-LIA Test Reagents	Preparation
Sample Diluent (<i>Green</i>)	Ready-to-use
Wash Solution (<i>Blue</i>)	1/5 in distilled water
Conjugate (<i>Red</i>)	Ready-to-use
Stop Solution (<i>White</i>)	Ready-to-use
Substrate (<i>Yellow</i>)	Ready-to-use
Drying Solution * (<i>Black</i>)	20% EtOH solution

* not included in the INNO-LIA Score kit. Should be prepared in the lab and has a 30-days stability (opened and unopened)


3. Start run and add samples

- Load samples in the sample carousel starting at position 1
Make sure the sample tubes are positioned straight in the carousel, with the barcodes facing outwards

REMARK: take into account a recommended dead volume of 1 mL and the required sample volume according to the assay IFU's


- Add the correct number of strips with the membrane side up and the strip label to the top, in the tray by using tweezers
- Load both reagent holders in the instrument and install the tubing's (wash and drying sol)
- Empty the tip waste container
- Load the tray with the strips in the instrument and close the holders
- Load the 1 mL tip rack and add 4 black 5 mL tips in the holder in the correct positions

- Empty the waste bottle




REMARK: After every action, tick the corresponding check box on the screen

- When running an **overnight protocol**, cover the tray with an adhesive sealer (the sealer added in the INNO-LIA kit can be used for this purpose) (see references)
- Place the sealer, above the markers covering appr. 80% of the troughs to prevent evaporation and leaving enough room on top of the tray to aspirate and dispense the wash solution




- *Incorrect positioning of the sealer, might lead to hampering of the instrument*
- Mark the checkbox to continue
- Push the "Scan" button to scan the sample barcodes (if option is activated)



REMARK: If the barcode scan option is deactivated, make sure each sample is loaded starting in carousel position 1 and in the correct order as required in the tray. If the same sample is required for multiple strips, add duplicates of this sample to each required carousel position.

- Press “Assay Start” to start the run



REMARK: Ensure that the instrument cover is correctly closed, otherwise the test will not proceed!



- After adding the samples, the instrument will prompt to check if the strips are still positioned correctly, facing upwards (remove the sample carousel if wanted)
- Press “Resume” to continue the run

4. Drying and scanning the strips

- When the test is done, the instrument will prompt you:
“Drying is about to start. Please remove the seal from the tray and ensure that the strips are positioned and oriented correctly, then press Resume to continue.”
- Remove the sealer, check that the strips are faced upwards and press resume
- Scanning of the strips will automatically start after one hour drying
- If running a LiRAS pre-generated worklist, the images will be automatically transferred to the interpretation software

III. Cleaning

- At the end of the assay, click on the cleaning icon to start the cleaning (section 5.2 in user manual)




- Load the cleaning solution in Wash and Drying bottles (4% Bleach or 4% SDS solution) and press “OK”
- After soaking, remove the cleaning solution, replace by deionized water and push “OK”
- Push “Done” after cleaning has finished



IV. Connection with LiRAS (IP-address)

When using the RoboBlot Device in combination with the LiRAS for Infectious Diseases software, follow these steps to obtain the IP address of the RoboBlot:

- Connect the instrument to the local network with the ethernet cable provided with the instrument
- On the instrument touch screen, press the button to open the options menu:



- Press and hold down on the Bee Robotics logo in the center of the screen. This activates Engineering Mode
- Navigate on the screen to **Engineering > General setting > Networking > Network service addresses:**



- Take note of the IP addresses and port number on the screen
- The device should have 2 IP-addresses. If the first address does not work, try the 2nd.

V. Calibration of the RoboBlot scanner

- When calibrating the RoboBlot scanner, use the 3 mm grey scale (see references) on the first position of a tray, with the label facing to the top of the tray
- Follow the instructions in LiRAS to calibrate the RoboBlot scanner

WARNING: The 3mm grey scale is only intended for calibration of RoboBlot devices, not for other scanning devices.

VI. References (FRE order codes)

- | | | | |
|---------------|--------------|------------------------------------|--------------------|
| - Trays: | 80742 | RoboBlot Tray 48 | 10 trays/pack |
| - 1 mL tips: | 80743 | 1000 µL Disposable tips (RoboBlot) | 10 trays (96 tips) |
| - 5 mL tips: | 80744 | 5000 µL Disposable tips (RoboBlot) | 10 trays (24 tips) |
| - Sealers: | 81061 | RoboBlot Sealers | 20 sealers/pack |
| - Grey scale: | 81100 | Grey scale package RoboBlot | 3 grey scales-3mm |

VII. Trademarks

- INNO-LIA is a trademark of Fujirebio Europe N.V., registered in US and other countries.
- RoboBlot is a trademark of Bee Robotics Ltd.